

A napari FF-SRM plugin for advanced microscopy image processing.


## 1. Installation process.

Review the system requirements before installing Anaconda Distribution. Select the installers for your computer's operating system:

<https://docs.anaconda.com/free/anaconda/install/index.html>


### 1.1 Installing Anaconda on Windows.

- 1) **Download** the Anaconda installer from <https://www.anaconda.com/> (fig. 1, left).
- 2) Go to the “Downloads” folder and double-click the **installer** to launch.
- 3) Click “**Next**” on the Setup window (fig. 1, right):




**Figure 1.** Installation windows for Anaconda Distribution. *Left:* Anaconda webpage. *Right:* Anaconda's setup window.

- 4) Read the licensing terms and click “**I Agree**” to continue (fig. 2, left).
- 5) It is recommended to install the program using the “**Just Me**” option, which will be limited to the current user account. However, if you wish the program available for all accounts on the computer, select “**All Users**”. Please note that selecting the latter option will require Windows Administrator permissions. Click “**Next**” (fig. 2, right):




**Figure 2.** Anaconda installation Setup. *Left:* License agreement. *Right:* Anaconda's installation type.

- 6) Select an installation **location** (Destination Folder) and click “**Next**” (fig. 3, left).
- 7) Choose whether to add Anaconda. It is recommended to register Anaconda software as your default Python because adding it to your PATH environment can interfere with other software. Click “**Install**” (fig. 3, right):



**Figure 3.** Anaconda installation Setup. *Left:* Installation location. *Right:* Anaconda's advanced installation options.

- 8) If you want to watch the packages Anaconda is installing, click “**Show details**”. Click “**Next**”.
- 9) After a successful installation you will see the “**Thanks for installing Anaconda**” dialog box.
- 10) Click the “**Finish**” button and verify the installation.



**Figure 4.** Anaconda installation Setup. *Left:* Installation details.

## 1.2 Installing Napari.

Napari is a fast, interactive and multi-dimensional image viewer for Python. To select which distribution to install and further details visit: <https://napari.org/stable/>.

For those familiar with Python, Napari can be installed on most macOS, Linux, and Windows systems using pip:

- 1) First install napari viewer: It is highly recommended to install napari into a clean virtual environment using an environment manager like [conda](#) or [venv](#). For example, with [conda](#):

```
conda create -y -n napari-env -c conda-forge python=3.9
conda activate napari-env
pip install "napari[all]"
```


- 2) Open Napari from command line:

Once installed, simply run  
napari

For further information visit: [https://napari.org/stable/tutorials/fundamentals/quick\\_start.html](https://napari.org/stable/tutorials/fundamentals/quick_start.html)


## 1.3 Installing napari-superres plugin.

- 1) When starting Napari, it will display the following window on the screen:



**Figure 5.** Napari's software window.

- 2) Select the “**Plugins**” Menu. Click the “**Install/Uninstall Plugins**” option (fig. 6).




**Figure 6.** Napari installation plugin menu.

### 3) Install napari-superres.


Simply search for the “napari-superres” plugin in the plugin install package interface and click install (fig. 7).

GitHub that contains the beta version of MSSR:  
<https://github.com/RoccoDAnt/napari-superres>




**Figure 7.** Search for the “napari-superres” (red rectangle) in the **Available Plugins** Napari menu to install.

### 4) Click “Install” and monitor the installation by clicking on **Show Status** option (fig. 8):



**Figure 8.** Napari-superres plugin installation process on Napari.

- 5) Once the plugin has been successfully installed, restart Napari and on the **“Plugins”** menu it will appear the napari-superres option and the following menu:



**Figure 9.** Napari-superres plugin installed on Napari.